



ALVERNO LABORATORIES

SPECIMEN COLLECTION GUIDE MOLECULAR

Specimen collection instructions and guidelines may change. We recommend checking the collection manual for the most up-to-date collection information at www.AlvernoLabs.com

Reviewed May 2023

Reach us at 800-937-5521

THIN PREP PAP



- ORDER# 144598
- ⇒ HPV
 - ⇒ GC/Chlamydia

ABBOTT ALINITY M MULTI-COLLECT SPECIMEN COLLEC-



- ORDER# 712363
- Intended for the collection and transportation of male and female swab and urine specimens for the detection of *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, & *Trichomonas vaginalis*

BBL CULTURE SWAB LIQUID STUART

ORDER# 193216



OR E-SWAB

ORDER# 132433

- ⇒ MRSA or MSSA by PCR



BD AFFIRM VAGINITIS PANEL

ORDER# 108401

- ⇒ Gardnerella
- ⇒ Trichomonas
- ⇒ Candida sp.*

*does NOT identify the species of Candida yeast



BD Max Vaginitis Panel

ORDER# 715839

- ⇒ Bacterial Vaginitis
- ⇒ Trichomonas
- ⇒ Candida species will be identified



E-SWAB

ORDER# 132433

- ⇒ Group B Strep by PCR
- ⇒ Vaginal/Rectal E-SWAB



STERILE SCREW CAP CUP (no additive)

ORDER# 153706

- ⇒ Respiratory Wash/Aspirate
 - * Respiratory 21 Pathogen Panel
 - * B. pertussis and parapertussis
 - * Mycoplasma pneumoniae
- ⇒ HSV 1 and 2
 - * Respiratory Wash / CSF
- ⇒ Varicella Zoster Virus (VZV)
 - * CSF
- ⇒ C. diff by PCR
 - * Stool
- ⇒ Enteric Pathogen Panels
 - * Stool (bacterial and parasitic)
- ⇒ Norovirus
 - * Stool
- ⇒ CMV Qualitative by PCR:
 - * CSF, Urine, Respiratory Wash
- ⇒ Meningitis/Encephalitis Pathogen Panel
 - * CSF



UTM (Universal Transport Media) ORDER# 123823

By PCR:

- ⇒ Varicella zoster virus (VZV)
- ⇒ Herpes virus 1&2
- ⇒ Mycoplasma pneumoniae
- ⇒ B pertussis and parapertussis
- ⇒ Respiratory 21 Pathogen Panel



UTM (Abbott) ORDER# 197838

- ⇒ SARS-CoV-2
- ⇒ Flu A / Flu B / RSV virus panel
- ⇒ Flu A / Flu B virus panel



LAVENDER TOP TUBE: WHOLE BLOOD

ORDER# 9603

- ⇒ Factor V Leiden Mutation by Real Time PCR
- ⇒ Factor II Mutation by Real Time PCR
- ⇒ CD4/CD8 (T-Cell Subsets)
- ⇒ Lymphocyte Subsets (T, B, and NK cells)
- ⇒ Lymphoma/Leukemia Immunophenotyping by FlowCytometry
- ⇒ Fetomaternal Bleed by FlowCytometry
- ⇒ HLA B27 Detection by FlowCytometry
- ⇒ Next Generation Sequencing, Myeloid Neoplasia Panel



UROVYSION

ORDER# 135884

- ⇒ For bladder cancer test collection use Urocyte Kit



PARA-PAK SAF FIXATIVE

ORDER# 11939

- ⇒ Enteric Parasite Panel

PARA-PAK C&S CARY-BLAIR MEDIA

ORDER# 23205

- ⇒ Gastrointestinal 22 Pathogen Panel
- ⇒ Enteric Bacterial Panel
- ⇒ Extended Enteric Bacterial Panel
- ⇒ Norovirus



LAVENDER TOP TUBE: PLASMA*

ORDER# 9603

- ⇒ CMV Quantitative by PCR
- ⇒ HCV Quantitative by PCR
- ⇒ HIV Quantitative by PCR
- ⇒ HBV Quantitative by PCR

*Submit plasma in a transfer tube



GREEN TOP TUBE: SODIUM HEPARIN

ORDER# 9604

- ⇒ Hematological FISH testing for peripheral blood or bone marrow aspirate

