

SPECIMEN COLLECTION GUIDE MOLECULAR

Specimen collection instructions and guidelines may change. We recommend checking the collection manual for the most up-to-date collection information at www.AlvernoLabs.com

Reviewed: May 2025 **Next Update: November 2025**

Reach us at 800-937-5521

THIN PREP PAP  ORDER# 144598 ⇒ HPV ⇒ GC/Chlamydia		ABBOTT ALINITY M MULTI-COLLECT SPECIMENT COLLEC-  ORDER# 712363 Intended for the collection and transportation of male and female swab and urine specimens for the detection of <i>Chlamydia trachomatis</i>, <i>Neisseria gonorrhoeae</i>, & <i>Trichomonas vaginalis</i>	
BBL CULTURE SWAB LIQUID STUART ORDER# 193216 OR E-SWAB ORDER# 132433 ⇒ MRSA or MSSA by PCR 		BD AFFIRM VAGINITIS PANEL ORDER# 108401 ⇒ Gardnerella ⇒ Trichomonas ⇒ Candida sp.* *does NOT identify the species of Candida yeast 	BD Max Vaginitis Panel ORDER# 715839 ⇒ Bacterial Vaginitis ⇒ Trichomonas ⇒ Candida species will be identified 
E-SWAB  ORDER# 132433 ⇒ Group B Strep by PCR ⇒ Vaginal/Rectal E-SWAB		STERILE SCREW CAP CUP/TUBE (no additive) ORDER# 153706 ⇒ Respiratory Wash/Aspirate * Respiratory 21 Pathogen Panel * B. pertussis and parapertussis * Mycoplasma pneumoniae ⇒ HSV 1 and 2 * Respiratory Wash / CSF ⇒ Varicella Zoster Virus (VZV) * CSF ⇒ C. diff by PCR * Stool ⇒ Enteric Pathogen Panels * Stool (bacterial and parasitic) ⇒ Norovirus * Stool ⇒ CMV Qualitative by PCR: * CSF, Urine, Respiratory Wash ⇒ Meningitis/Encephalitis Pathogen Panel * CSF 	
UTM (Universal Transport Media) ORDER# 123823 By PCR: ⇒ Varicella zoster virus (VZV) ⇒ Herpes virus 1&2 ⇒ Mycoplasma pneumoniae ⇒ B pertussis and parapertussis ⇒ Respiratory 21 Pathogen Panel 		PARA-PAK SAF FIXATIVE ORDER# 11939 ⇒ Enteric Parasite Panel PARA-PAK C&S CARY-BLAIR MEDIA ORDER# 23205 ⇒ Gastrointestinal 22 Pathogen Panel ⇒ Enteric Bacterial Panel ⇒ Extended Enteric Bacterial Panel ⇒ Norovirus 	
UTM (Abbott) ORDER# 197838 ⇒ SARS-CoV-2 ⇒ Flu A / Flu B / RSV virus panel ⇒ Flu A / Flu B virus panel 			
LAVENDER TOP TUBE: WHOLE BLOOD ORDER# 9603 ⇒ Factor V Leiden Mutation by Real Time PCR ⇒ Factor II Mutation by Real Time PCR ⇒ CD4/CD8 (T-Cell Subsets) ⇒ Lymphocyte Subsets (T, B, and NK cells) ⇒ Lymphoma/Leukemia Immunophenotyping by FlowCytometry ⇒ Fetomaternal Bleed by FlowCytometry ⇒ HLA B27 Detection by FlowCytometry ⇒ Next Generation Sequencing, Myeloid Neoplasia Panel 			
LAVENDER TOP TUBE: PLASMA* ORDER# 9603 ⇒ CMV Quantitative by PCR ⇒ HCV Quantitative by PCR ⇒ HIV Quantitative by PCR ⇒ HBV Quantitative by PCR *Submit plasma in a transfer tube 		GREEN TOP TUBE: SODIUM HEPARIN ORDER# 9604 ⇒ Hematological FISH testing for peripheral blood or bone marrow aspirate 	