














SPECIMEN COLLECTION GUIDE MOLECULAR

<p>THIN PREP PAP ORDER# 144598</p>	<p>Abbott Multi-Collect Specimen Collection Kit</p>	
 <ul style="list-style-type: none"> ⇒ HPV ⇒ GC/Chlamydia 		<p>ORDER# 176526 Intended for the collection and transportation of male and female, swab and urine specimens for the detection of Chlamydia trachomatis and Neisseria gonorrhoeae</p>
<p>BBL Culture Swab Liquid Stuart OR</p>  <p>ORDER# 193216</p> <p>E-SWAB ORDER# 132433</p> <ul style="list-style-type: none"> ⇒ MRSA/MSSA by PCR ⇒ MRSA by PCR 	<p>BD Affirm Vaginitis Panel</p> <p>ORDER# 108401</p>  <ul style="list-style-type: none"> ⇒ Gardnerella ⇒ Trichomonas ⇒ Candida sp. 	<p>BD Max Vaginitis Panel</p> <p>ORDER# 181279</p>  <ul style="list-style-type: none"> ⇒ Bacterial Vaginitis ⇒ Trichomonas ⇒ Candida sp., C. glabrata, C. Krusei
<p>E-SWAB</p> <p>ORDER# 132433</p> <ul style="list-style-type: none"> ⇒ Group B Strep by PCR <p>Vaginal/Rectal eSwab</p> 	<p>STERILE SCREW CAP CUP (no additive)</p> <p>ORDER# 153706</p> <ul style="list-style-type: none"> ⇒ Respiratory Wash/Aspirate <ul style="list-style-type: none"> * Flu A / Flu B / RSV panel * Flu A / Flu B panel * Respiratory 21 Pathogen Panel * B. pertussis and parapertussis * Mycoplasma pneumoniae ⇒ HSV 1 and 2 <ul style="list-style-type: none"> * Respiratory Wash / CSF ⇒ Varicella Zoster Virus (VZV) <ul style="list-style-type: none"> * CSF ⇒ C. diff by PCR <ul style="list-style-type: none"> * Stool ⇒ Enteric Pathogen Panels <ul style="list-style-type: none"> * Stool (bacterial and parasitic) ⇒ Norovirus <ul style="list-style-type: none"> * Stool ⇒ CMV Qualitative by PCR: <ul style="list-style-type: none"> * CSF, Urine, Respiratory Wash 	
<p>BD UTM (Universal Transport Media)</p> <p>ORDER# 111005</p> <p>By PCR:</p> <ul style="list-style-type: none"> ⇒ Varicella virus ⇒ Herpes virus 1&2 ⇒ Flu A / Flu B / RSV virus panel ⇒ Flu A / Flu B virus panel ⇒ Mycoplasma pneumoniae ⇒ B pertussis and parapertussis ⇒ Respiratory 21 Pathogen Panel 	<p>UroVysion ORDER# 135884</p> <ul style="list-style-type: none"> ⇒ For bladder cancer test collection use Urocyte Kit 	
<p>PURPLE TOP TUBE - WHOLE BLOOD</p> <p>ORDER# 9603</p> <ul style="list-style-type: none"> ⇒ Factor V Leiden Mutation by Real Time PCR ⇒ Factor II Mutation by Real Time PCR ⇒ Cystic Fibrosis screen: original tube only ⇒ CD4/CD8 (T-Cell Subsets) ⇒ Lymphocyte Subsets (T, B, and NK cells) ⇒ Lymphoma/Leukemia Immunophenotyping by FlowCytometry ⇒ Fetomaternal Bleed by FlowCytometry ⇒ HLA B27 Detection by FlowCytometry 	<p>Para-Pak SAF Fixative ORDER# 11939</p> <ul style="list-style-type: none"> ⇒ Enteric Parasite Panel <p>Para-Pak C&S Cary-Blair Media ORDER# 23205</p> <ul style="list-style-type: none"> ⇒ Gastrointestinal 22 Pathogen Panel ⇒ Enteric Bacterial Panel ⇒ Extended Enteric Bacterial Panel ⇒ Norovirus  	
<p>PURPLE TOP TUBE - PLASMA</p> <p>ORDER# 9603</p> <ul style="list-style-type: none"> ⇒ CMV Qualitative PCR ⇒ CMV Quantitative PCR ⇒ HCV Quantitative by PCR ⇒ HIV Quantitative by PCR ⇒ HBV Quantitative By PCR 	<p>GREEN TOP TUBE SODIUM HEPARIN</p> <p>ORDER# 9604</p> <ul style="list-style-type: none"> ⇒ Hematological FISH testing for peripheral blood or bone marrow aspirate 